

DR. F. H. H. BRILL · C/O DR. BRILL + PARTNER GMBH · STIEGSTÜCK 34 · DE-22339 HAMBURG

ROPIMEX R. OPEL GmbH
Bildstocker Straße 12
DE - 66538 Neunkirchen

Hamburg, 27 January 2022

Expert opinion

Bactericidal Activity of **Bacoban WBUS** in the quantitative suspension test according to DIN EN 1656:2019 (Phase 2, Step 1)

The disinfectant **Bacoban WBUS** was tested and evaluated according to DIN EN 1656:2019 "Chemical disinfectants and antiseptics - Quantitative suspension test for the evaluation of bactericidal activity of chemical disinfectants and antiseptics used in veterinary area - Test method and requirements (phase 2, step 1)".

According to the test report no. L21/01041.5 dated 27/01/2022 of Dr. Brill + Partner GmbH the preparation showed bactericidal activity under clean conditions at a test temperature of $10^{\circ}\text{C} \pm 1^{\circ}\text{C}$.

Bacoban WBUS complies with the requirements of DIN EN 1656:2019 (phase 2, step 1) with the following concentration-time relationship:

Bactericidal:	clean conditions	0.5 %	30 minutes
----------------------	-------------------------	--------------	-------------------



Dr. Florian H. H. Brill

Test report no L21/01041.5

Quantitative suspension test for the evaluation of bactericidal activity of **Bacoban WBUS**
in the veterinary area (DIN EN 1656:2019; Phase 2, Step 1*)

In accordance with your order, we tested the preparation **Bacoban WBUS** for its activity in the quantitative suspension test according to DIN EN 1656:2019* under clean conditions.

1 General Information and Material

1.1 Client

Client: ROPILEX R. OPEL GmbH, Mrs Jennifer Sahl, Bildstocker Straße 12,
DE - 66538 Neunkirchen, Germany
Date of order: 23/12/2021
Confirmation no.: 226415

1.2 Identification of Test Laboratory

Location: Dr. Brill + Partner GmbH · Institute for Hygiene and Microbiology,
Stiegstück 34, DE-22339 Hamburg, Germany
Study manager: Dipl.-Ing. Dr. rer. nat. Andreas Kampe
Scientific assistant: Dipl.-Biol. Henrik Gabriel
Laboratory technicians: Elahe Saroukhani

1.3 Table of Contents

General Information and Material	1
Methods	3
Results	3
Tables	4
List of Abbreviations	8

1.4 Identification of Sample

Name of product: **Bacoban WBUS**
Batch no.: 20210920_Bacoban WBUS_imi

* Test procedure accredited according to DIN EN ISO/IEC 17025. Test report issued by Dr. Brill + Partner GmbH, Stiegstück 34, DE - 22339 Hamburg, Phone +49 40 557631-0, Telefax +49 40 557631-11, www.brillhygiene.com. No copying or transmission, in whole or in part, of this test report without the explicit prior written permission. The test results exclusively apply to the tested samples. Information on measurement uncertainty and Version history on request. © Dr. Brill + Partner GmbH 2022



Internal no.:	21/01202
Manufacturer:	ROPIMEX R. OPEL GmbH, DE - 66538 Neunkirchen, Germany
Date of delivery:	21/09/2021
Storage conditions:	room temperature and darkness
Appearance of product:	clear liquid
Odour:	characteristic
Product type:	surface disinfectant
Recommended diluent:	Tap water
Diluent used:	water of standardised hardness (WSH, pH 7.0)
pH value, concentrate:	4.3
pH value, 1.0 % (measured in WSH):	6.2
pH value, 0.5 % (measured in WSH):	6.6
pH value, 0.1 % (measured in WSH):	7.0
Active agents (Manufacturer's data):	33.31 g benzalkoniumchloride

1.5 Test Conditions

Test period:	04/01/ - 07/01/2022
Lab task no.:	L21/01041.5
Product test concentrations:	0.1 + 0.5 + 1.0 %
Exposure time:	30 minutes
Test temperature:	10°C ± 1°C
Incubation temperature:	36°C ± 1°C
Organic load:	clean conditions (3.0 g/L bovine albumin)
Neutraliser:	60 g/L polysorbate 80, 60 g/L saponine, 8 g/L lecithin, 1 g/L cysteine, 2.5 g/L SDS (TLSC-SDS)
Test organisms:	<i>Staphylococcus aureus</i> ATCC 6538 <i>Enterococcus hirae</i> ATCC 10541 <i>Proteus hauseri (vulgaris)</i> ATCC 13315 <i>Pseudomonas aeruginosa</i> ATCC 15442

* Test procedure accredited according to DIN EN ISO/IEC 17025. Test report issued by Dr. Brill + Partner GmbH, Stiegstück 34, DE - 22339 Hamburg, Phone +49 40 557631-0, Telefax +49 40 557631-11, www.brillhygiene.com. No copying or transmission, in whole or in part, of this test report without the explicit prior written permission. The test results exclusively apply to the tested samples. Information on measurement uncertainty and Version history on request. © Dr. Brill + Partner GmbH 2022



2 Methods

The tests were carried out according to DIN EN 1656:2019 "Chemical disinfectants and antiseptics - Quantitative suspension test for the evaluation of bactericidal activity of chemical disinfectants and antiseptics used in veterinary area - Test method and requirements (phase 2, step 1)".


3 Results

The test results based on DIN EN 1656:2019 are summarised in tables 1.

The test bacteria were sufficiently (RF >5) inactivated with the following concentration-time relationship:

Bactericidal: clean conditions 0.5 % 30 minutes

Hamburg, 27/01/2022



Dipl.-Ing. Dr. rer. nat. Andreas Kampe
Study Manager



Dipl.-Biol. Henrik Gabriel
Head of Laboratory



* Test procedure accredited according to DIN EN ISO/IEC 17025. Test report issued by Dr. Brill + Partner GmbH, Stiegstück 34, DE - 22339 Hamburg, Phone +49 40 557631-0, Telefax +49 40 557631-11, www.brillhygiene.com. No copying or transmission, in whole or in part, of this test report without the explicit prior written permission. The test results exclusively apply to the tested samples. Information on measurement uncertainty and Version history on request. © Dr. Brill + Partner GmbH 2022



Table 1.1: Validation, Controls and Evaluation

Product name: **Bacoban WBUS** Batch: 20210920_Bacoban
WBUS_imi
Test organism: *Staphylococcus aureus* Temperature: 10°C ± 1°C
Organic load: clean conditions Neutraliser: TLSC-SDS

Suspension for validation (N _{v0})			Control of test conditions (A)						
			5 minutes		30 minutes		n.t.		
Microbial count	∑		Microbial count	∑	Microbial count	∑	Microbial count	∑	
V _{c1}	113	116,5	V _{c1}		V _{c1}	117	120,5	V _{c1}	
V _{c2}	120		V _{c2}		V _{c2}	124		V _{c2}	
30 ≤ ∑ of N _{v0} ≤ 160		Yes	∑ of A(5') is ≥ 0,5 x ∑ of N _{v0} ?		∑ of A(30') is ≥ 0,5 x ∑ of N _{v0} ?		Yes	∑ of A(60') is ≥ 0,5 x ∑ of N _{v0} ?	
Control of neutraliser (B)			Validation of method (C) at highest product concentration: 1,0 %						
			5 minutes		30 minutes		n.t.		
Microbial count	∑		Microbial count	∑	Microbial count	∑	Microbial count	∑	
V _{c1}	109	113	V _{c1}		V _{c1}	101	103	V _{c1}	
V _{c2}	117		V _{c2}		V _{c2}	105		V _{c2}	
∑ of B is ≥ 0,5 x ∑ of N _{v0} ?		Yes	∑ of C(5') is ≥ 0,5 x ∑ of N _{v0} ?		∑ of C(30') is ≥ 0,5 x ∑ of N _{v0} ?		Yes	∑ of C(60') is ≥ 0,5 x ∑ of N _{v0} ?	

Test suspension (N and N ₀)	N	Microbial count of plates		V _{c1}	V _{c2}	∑ _{wm} / lg N	N ₀ =N/10; lg N ₀	7,17 ≤ N ₀ ≤ 7,70 ?
	1,00E-06	>330	>330	>330	>330	4,40E+08	7,64	Yes
	1,00E-07	44	44	44	44	8,64		
Product concentration [%]	Exposure time [min]	Microbial count of plates		V _{c1}	V _{c2}	N _a = ∑ x 10	lg N _a	lg R (lg N ₀ = 7,64)
0,10	30	0	0	<14	<14	<140	<2,15	≥ 5,49
0,50	30	0	0	<14	<14	<140	<2,15	≥ 5,49
1,00	30	0	0	<14	<14	<140	<2,15	≥ 5,49

* Test procedure accredited according to DIN EN ISO/IEC 17025. Test report issued by Dr. Brill + Partner GmbH, Stiegstück 34, DE - 22339 Hamburg, Phone +49 40 557631-0, Telefax +49 40 557631-11, www.brillhygiene.com. No copying or transmission, in whole or in part, of this test report without the explicit prior written permission. The test results exclusively apply to the tested samples. Information on measurement uncertainty and Version history on request. © Dr. Brill + Partner GmbH 2022



Table 1.2: Validation, Controls and Evaluation

Product name: **Bacoban WBUS** Batch: 20210920_Bacoban
WBUS_imi
Test organism: *Enterococcus hirae* Temperature: 10°C ± 1°C
Organic load: clean conditions Neutraliser: TLSC-SDS

Suspension for validation (N _{v0})				Control of test conditions (A)											
				5 minutes				30 minutes				n.t.			
Microbial count		x̄		Microbial count		x̄		Microbial count		x̄		Microbial count		x̄	
V _{c1}	84		86,5	V _{c1}				V _{c1}	80		81,5	V _{c1}			
V _{c2}	89			V _{c2}				V _{c2}	83			V _{c2}			
30 ≤ x̄ of N _{v0} ≤ 160				Yes				x̄ of A(5') is ≥ 0,5 x x̄ of N _{v0} ?				Yes			
				x̄ of A(30') is ≥ 0,5 x x̄ of N _{v0} ?								x̄ of A(60') is ≥ 0,5 x x̄ of N _{v0} ?			
Control of neutraliser (B)				Validation of method (C) at highest product concentration: 1,0 %											
				5 minutes				30 minutes				n.t.			
Microbial count		x̄		Microbial count		x̄		Microbial count		x̄		Microbial count		x̄	
V _{c1}	79		88	V _{c1}				V _{c1}	77		79,5	V _{c1}			
V _{c2}	97			V _{c2}				V _{c2}	82			V _{c2}			
x̄ of B is ≥ 0,5 x x̄ of N _{v0} ?				Yes				x̄ of C(5') is ≥ 0,5 x x̄ of N _{v0} ?				Yes			
				x̄ of C(30') is ≥ 0,5 x x̄ of N _{v0} ?								x̄ of C(60') is ≥ 0,5 x x̄ of N _{v0} ?			
Test suspension (N and N ₀)		N		Microbial count of plates		V _{c1}		V _{c2}		x̄ _{wm} / lg N		N ₀ =N/10; lg N ₀		7,17 ≤ N ₀ ≤ 7,70 ?	
		1,00E-06		>330		>330		>330		3,30E+08		7,52		Yes	
		1,00E-07		29		37		29		8,52					
Product concentration [%]		Exposure time [min]		Microbial count of plates		V _{c1}		V _{c2}		N _a = x̄ x 10		lg N _a		lg R (lg N ₀ = 7,52)	
0,10		30		0		<14		<14		<140		<2,15		≥ 5,37	
0,50		30		0		<14		<14		<140		<2,15		≥ 5,37	
1,00		30		0		<14		<14		<140		<2,15		≥ 5,37	

* Test procedure accredited according to DIN EN ISO/IEC 17025. Test report issued by Dr. Brill + Partner GmbH, Stiegstück 34, DE - 22339 Hamburg, Phone +49 40 557631-0, Telefax +49 40 557631-11, www.brillhygiene.com. No copying or transmission, in whole or in part, of this test report without the explicit prior written permission. The test results exclusively apply to the tested samples. Information on measurement uncertainty and Version history on request. © Dr. Brill + Partner GmbH 2022



Table 1.3: Validation, Controls and Evaluation

Product name: **Bacoban WBUS** Batch: 20210920_Bacoban
WBUS_imi
Test organism: *Proteus hauseri* Temperature: 10°C ± 1°C
Organic load: clean conditions Neutraliser: TLSC-SDS

Suspension for validation (N _{v0})			Control of test conditions (A)											
			5 minutes				30 minutes				n.t.			
Microbial count		̄	Microbial count		̄	Microbial count		̄	Microbial count		̄	Microbial count		̄
V _{c1}	103	110	V _{c1}		92	V _{c1}	84	92	V _{c1}			V _{c1}		
V _{c2}	117		V _{c2}			V _{c2}	100		V _{c2}					
30 ≤ ̄ of N _{v0} ≤ 160		Yes	̄ of A(5') is ≥ 0,5 x ̄ of N _{v0} ?				̄ of A(30') is ≥ 0,5 x ̄ of N _{v0} ?				Yes ̄ of A(60') is ≥ 0,5 x ̄ of N _{v0} ?			
Control of neutraliser (B)			Validation of method (C) at highest product concentration: 1,0 %											
			5 minutes				30 minutes				n.t.			
Microbial count		̄	Microbial count		̄	Microbial count		̄	Microbial count		̄	Microbial count		̄
V _{c1}	95	95,5	V _{c1}		93	V _{c1}	89	93	V _{c1}			V _{c1}		
V _{c2}	96		V _{c2}			V _{c2}	97		V _{c2}					
̄ of B is ≥ 0,5 x ̄ of N _{v0} ?		Yes	̄ of C(5') is ≥ 0,5 x ̄ of N _{v0} ?				̄ of C(30') is ≥ 0,5 x ̄ of N _{v0} ?				Yes ̄ of C(60') is ≥ 0,5 x ̄ of N _{v0} ?			
Test suspension (N and N ₀)		N	Microbial count of plates		V _{c1}	V _{c2}	̄ _{wm} / lg N	N ₀ =N/10; lg N ₀	7,17 ≤ N ₀ ≤ 7,70 ?					
		1,00E-06	>330	>330	>330	>330	4,70E+08	7,67	Yes					
		1,00E-07	46	48	46	48	8,67							
Product concentration [%]		Exposure time [min]	Microbial count of plates		V _{c1}	V _{c2}	N _a = ̄ x 10	lg N _a	lg R (lg N ₀ = 7,67)					
0,10		30	0	0	<14	<14	<140	<2,15	≥ 5,52					
0,50		30	0	0	<14	<14	<140	<2,15	≥ 5,52					
1,00		30	0	0	<14	<14	<140	<2,15	≥ 5,52					

* Test procedure accredited according to DIN EN ISO/IEC 17025. Test report issued by Dr. Brill + Partner GmbH, Stiegstück 34, DE - 22339 Hamburg, Phone +49 40 557631-0, Telefax +49 40 557631-11, www.brillhygiene.com. No copying or transmission, in whole or in part, of this test report without the explicit prior written permission. The test results exclusively apply to the tested samples. Information on measurement uncertainty and Version history on request. © Dr. Brill + Partner GmbH 2022



Table 1.4: Validation, Controls and Evaluation

Product name: **Bacoban WBUS** Batch: 20210920_Bacoban
Test organism: *Pseudomonas aeruginosa* WBUS_imi
Organic load: clean conditions Temperature: 10°C ± 1°C
Neutraliser: TLSC-SDS

Suspension for validation (N _{v0})			Control of test conditions (A)						
			5 minutes		30 minutes		n.t.		
Microbial count	x̄		Microbial count	x̄	Microbial count	x̄	Microbial count	x̄	
V _{c1}	85	85	V _{c1}		V _{c1}	83	87,5	V _{c1}	
V _{c2}	85		V _{c2}		V _{c2}	92		V _{c2}	
30 ≤ x̄ of N _{v0} ≤ 160		Yes	x̄ of A(5') is ≥ 0,5 x x̄ of N _{v0} ?		x̄ of A(30') is ≥ 0,5 x x̄ of N _{v0} ?		Yes	x̄ of A(60') is ≥ 0,5 x x̄ of N _{v0} ?	
Control of neutraliser (B)			Validation of method (C) at highest product concentration: 1,0 %						
			5 minutes		30 minutes		n.t.		
Microbial count	x̄		Microbial count	x̄	Microbial count	x̄	Microbial count	x̄	
V _{c1}	84	89,5	V _{c1}		V _{c1}	92	92,5	V _{c1}	
V _{c2}	95		V _{c2}		V _{c2}	93		V _{c2}	
x̄ of B is ≥ 0,5 x x̄ of N _{v0} ?		Yes	x̄ of C(5') is ≥ 0,5 x x̄ of N _{v0} ?		x̄ of C(30') is ≥ 0,5 x x̄ of N _{v0} ?		Yes	x̄ of C(60') is ≥ 0,5 x x̄ of N _{v0} ?	

Test suspension (N and N ₀)	N	Microbial count of plates		V _{c1}	V _{c2}	x̄ _{wm} / lg N	N ₀ =N/10; lg N ₀	7,17 ≤ N ₀ ≤ 7,70 ?
	1,00E-06	>330	>330	>330	>330	4,60E+08	7,66	Yes
	1,00E-07	43	49	43	49	8,66		
Product concentration [%]	Exposure time [min]	Microbial count of plates		V _{c1}	V _{c2}	N _a = x̄ x 10	lg N _a	lg R (lg N ₀ = 7,66)
0,10	30	>330	>330	>330	>330	>3300	>3,52	≤ 4,14
0,50	30	20	21	20	21	205	2,31	5,35
1,00	30	0	0	<14	<14	<140	<2,15	≥ 5,51

* Test procedure accredited according to DIN EN ISO/IEC 17025. Test report issued by Dr. Brill + Partner GmbH, Stiegstück 34, DE - 22339 Hamburg, Phone +49 40 557631-0, Telefax +49 40 557631-11, www.brillhygiene.com. No copying or transmission, in whole or in part, of this test report without the explicit prior written permission. The test results exclusively apply to the tested samples. Information on measurement uncertainty and Version history on request. © Dr. Brill + Partner GmbH 2022



4 List of Abbreviations

A	=	control of test conditions
B	=	control of neutraliser
C	=	validation of method at highest product concentration
N	=	test suspension
N _{vo}	=	suspension for validation
n.t.	=	not tested
N ₀	=	microbial count of test suspension N / 10 (microbial count at time index 0)
R	=	germ reduction in log ₁₀ -steps
V _c	=	viable microbial count per ml
\bar{x}	=	weighted mean of N

* Test procedure accredited according to DIN EN ISO/IEC 17025. Test report issued by Dr. Brill + Partner GmbH, Stiegstück 34, DE - 22339 Hamburg, Phone +49 40 557631-0, Telefax +49 40 557631-11, www.brillhygiene.com. No copying or transmission, in whole or in part, of this test report without the explicit prior written permission. The test results exclusively apply to the tested samples. Information on measurement uncertainty and Version history on request. © Dr. Brill + Partner GmbH 2022