



Shelf-life of Bacoban® products

Product	Art. No.	Unopened	Opened
WB concentrate	BACWB1D BACWB5	24 months	24 months (appropriate handling and storage provided)
WB, diluted with tap water	BACWB1 BACWB5	8 hours	
WB-DL concentrate	BACWBDL	12 months	6 months
WB-DL diluted with demineralized water (<math><0,11\mu\text{S}/\text{cm}</math> conductivity, <math>0^\circ\text{dh}< <math="" math>,="">\leq 100\text{ CFU}/\text{ml}</math>)</math>0^\circ\text{dh}<>	BACD500DI X-BAC500DL1 BAC1DL3	12-36 months	Spray bottle: 12 months Canister with screw cap: 36 months (appropriate handling and storage provided)
Wet wipes in flowpack	BACDLTUE	36 months	1 month
Alcohol based	BAC500 spray BAC5000 canister	12-24 months	Spray bottle: 12 months Canister with screw cap: 36 months (appropriate handling and storage provided)